


# The Welch Allyn Green Series™ 900 Procedure Light



6 LEDs



Touchless On/Off  
with Dimming Control

Disposable Sheath  
Available—Use #52630

Green Series™ Medical Exam Lights



WelchAllyn®

Advancing Frontline Care™

# Whiter, Brighter, and Greener

The GS 900 Procedure Light features:

- Six LEDs, combining outstanding durability, extraordinary maneuverability, unparalleled spot quality and multiple mounting options
- Minimal heat output and high color temperature
- Less energy consumption
- No bulb replacement



## IDEAL FOR USE IN THE FOLLOWING ENVIRONMENTS:

- AMBULATORY SURGERY CENTERS
- HOSPITALS—ER AND ICU
- OB/GYN
- PLASTIC SURGERY
- MINOR SURGERY OUTPATIENT FACILITIES

## MOUNTING OPTIONS AVAILABLE:



|                           |                                    |                                                                                                     |
|---------------------------|------------------------------------|-----------------------------------------------------------------------------------------------------|
| Light Specifications      | # of LEDs                          | 6 White LEDs                                                                                        |
|                           | Colour Temperature                 | 5,500K                                                                                              |
|                           | Light Intensity @ Working Distance | 80,000 Lux @ 61cm (24in)<br>65,000 Lux @ 100cm (40in)                                               |
|                           | Spot Size @ Working Distance       | 28cm @ 61cm (24in)<br>23cm @ 100cm (40in)                                                           |
| Physical Dimensions       | Reach Mobile Stand                 | 95cm (37,5in)                                                                                       |
|                           | Reach Wall                         | 175cm (69in)                                                                                        |
|                           | Reach Ceiling                      | 175cm (69in)                                                                                        |
|                           | Ceiling up/down                    | 95cm (37,4in)                                                                                       |
|                           | Luminaire (Height/Width/Depth)     | 34,9 / 45,8 / 10,2 cm<br>(13,74 / 18,03 / 4,02in)                                                   |
|                           | Control Box (Height/Width/Depth)   | n/a                                                                                                 |
|                           | Mobile Stand (Height / Diameter)   | Mobile Stand Base:<br>Length: 55,8cm (21,97in)<br>Width: 56,5cm (22,24in)<br>Height: 16,5cm (6,5in) |
|                           | Total Weight Table/Wall mount      | 8,12 kg (17,90 lbs)                                                                                 |
| Total Weight Mobile Stand | 20,34 kg (44,85 lbs)               |                                                                                                     |
| Electrical                | Power Supply                       | Input: 100-240 V~ 50Hz-60Hz,<br>1,00A max                                                           |
|                           | Power Cord                         | #18 AWG, 3-wire grounded cord set,<br>8' (2,5 m) long                                               |
|                           | Control Functions                  | Touch less On/Off with Dimming                                                                      |
| Misc Data                 | Mounting Options                   | Mobile Stand, Wall mount, Ceiling mount                                                             |
|                           | Kind of Extension                  | Flexible arm system with obedient articulation                                                      |
|                           | Accessory Options                  | Disposable sheath: P/N 52630<br>Length: 22,86 cm (9in)                                              |
|                           | Lifetime                           | 50,000 hours                                                                                        |
|                           | Warranty                           | 5-year limited warranty                                                                             |

To place an order, contact your appropriate Customer Service Department.



Welch Allyn Corporate Headquarters  
4341 State Street Road, P.O. Box 220  
Skaneateles Falls, NY 13153-0220, U.S.A.  
Telephone: +1.315.685.4602  
+1.800.535.6663  
Fax: +1.315.685.0054

Welch Allyn European Customer Service: +353 46 90 67790  
Australia: +61 29 638 3000  
China: +86 21 6327 9631  
France: +33 160 09 33 66  
France (Customer Service): +33 155 69 58 49  
Germany: +49 6950 985 132  
Spain: +34 912 990 376  
Japan: +81 3 3219 0071  
India: +91 44 4353 5724

Latin America (Miami office): +1 305 669 9003  
Malaysia: +603 78753341  
Netherlands: +31 202 061 360  
Singapore: +65 6419 8100  
Sweden: +46 042 33 02 20  
United Kingdom: +44 207 365 6780

[www.welchallyn.com/greenserries](http://www.welchallyn.com/greenserries)

**WelchAllyn**

Advancing Frontline Care™

