

Primo™

A simple-to-use fully automatic therapy mattress that responds to the patient



Features and benefits

- Specifically designed for enhanced prevention and treatment of pressure ulcers
- Hill-Rom I-mmersion™ Therapy automatically redistributes pressure to a larger surface area, significantly reducing compression of vulnerable tissues
- Internal I-mmersion sensor automatically adjusts cell pressures to ensure the lowest possible interface pressures are achieved at all times
- Compact and lightweight air supply unit is extremely quiet, vibration free and very user friendly
- One touch P-MAX feature provides fast inflation to aid patient repositioning and transfer

Primo

Dynamic low pressure therapy system

The Primo mattress replacement system from Hill-Rom helps to prevent and treat pressure ulcers in low to high risk adult patients by providing continuous low pressure redistribution.

The precise level of pressure redistribution is continuously and automatically regulated by the internal immersion sensor which adjusts internal cell pressures according to each patient's size, weight, morphology and position.

Twenty-four-hour pressure ulcer management can be achieved with the addition of the Reflex seat cushion which connects to the same air supply unit, yet can be sealed for greater patient independence whilst out of bed.



Technical Specifications

Mattress Replacement

Length	200 cm
Width	90 cm
Width (narrow version)	77cm
Height	16 cm
Weight	8.4 kg

Control Unit

Length	29 cm
Width	15 cm
Height	12 cm
Weight	3.5 kg

Safe Working Load	250 kg
-------------------	--------

Recommended therapeutic weight	30 to 150 kg
--------------------------------	--------------

Cover material

Polyurethane coating on polyester material, 360° zipper, welded seams, low friction, two-way stretch, moisture vapour permeable, bacteriostatic, fungistatic and antimicrobial. Can be wiped and washed. Chlorine-based cleaning solutions may be used ($\leq 1,000$ ppm).

Hill-Rom is a leading global medical technology company with approximately 10,000 employees worldwide. We partner with health care providers in more than 100 countries by focusing on patient care solutions that improve clinical and economic outcomes in five core areas: **Advancing Mobility, Wound Care and Prevention, Patient Monitoring and Diagnostics, Surgical Safety and Efficiency, and Respiratory Health.** Around the world the people, programs, and products of Hill-Rom work towards one mission: **Enhancing outcomes for patients and their caregivers.**

Not all products/options are available in all countries. For further information about our products and services, please contact your local Hill-Rom representative or visit our webpage:

www.hill-rom.com

The Primo surface is intended to help prevent and treat pressure ulcers Stage I, II and III in low to high risk patients, with a minimum weight of 30 kg and a minimum height of 1.46 m, in a range of healthcare facilities.

Class IIa

Notified Body: LNE/G-Med CE0459

Manufacturer: Hill-Rom S.A.S. – BP 14 ZI du Talhouët 56330 Pluvigner France

This medical device is a regulated health product which, pursuant to such regulation bears a CE mark. Hill-Rom recommends that you carefully read the detailed instructions for safe and proper use included in the documents accompanying the medical devices. The personnel of healthcare establishments are responsible for the proper use and maintenance of these medical devices. Hill-Rom reserves the right to make changes without notice in design, specifications and models. The only warranty Hill-Rom makes is the express written warranty extended on the sale or rental of its products.

©2017 Hill-Rom Services, Inc. ALL RIGHTS RESERVED.

Doc. No: 5EN119304-03, 01 May 2017

Enhancing outcomes for patients and their caregivers:

Hill-Rom